MACHEREY-NAGEL Chromatography application note

MN application note 07/2020

Solid phase extraction of per- and polyfluoroalkyl substances (PFAS) from contaminated soils

MACHEREY-NAGEL application department · Dr. H. R. Wollseifen, T. Kretschmer, S. Schneider

Abstract

This application note describes the determination of per- and polyfluoroalkyl substances (PFAS) from contaminated soils. It demonstrates the extraction of PFAS from soil samples using CHROMABOND[®] PFAS column, a special SPE combination phase, for the methodology described in DIN 38407-42. The eluates are finally analyzed by HPLC-MS/MS.

Introduction

Per- and polyfluoroalkyl substances (PFAS) have been manufactured since the 1940s and have been used for various applications due to their unique chemical properties. PFAS are used as additives for example in:

- Fire-fighting foam
- Fiber coating
- Textile coating, e.g. seat covers, carpets, outdoor clothing
- Cookware
- Paper finishing
- Food packaging, e.g. pizza cartons, paper cups
- Building material, e.g. water resistant lacquer.

This broad use, appearance and their persistency leads to the fact that PFAS are now abundant in the environment. From fire-fighting foams, textilies and food-packaging PFAS blaze their trail into soil and from there into ground water. During this distribution, degradation processes can also take place leading to several PFAS by-products. From the ground water, PFAS and their degradation products are being further transported to other locations in the environment. Unfortunately, many of the PFAS are toxic. Therefore, monitoring these substances is important.



Figure 1: General structure of various per- and polyfluoroalkyl compounds (PFAS). In a previous application note [1], we presented a solution for the analysis of PFAS from water. CHROMABOND[®] PFAS solid phase extraction columns showed excellent recovery rates and reproducibility for 30 PFAS. This present application note describes the use of CHROMABOND[®] PFAS for the enrichment of several PFAS from different types of soil according to DIN 38414-14 [2]. In addition to

the 10 analytes mentioned in DIN 38414-14, the note even describes the enrichment of a total of 40 PFAS. The extracts are analyzed by HPLC-MS/MS.

Sample pretreatment for solid phase extraction (SPE)

Weigh out 2.5 g of homogenized sample (dried) into a 50 mL centrifuge tube

- Add 62.5 μ L of standard solution (β = 0.2 μ g/mL for each compound in methanol) for determining recovery rate
- Add 25 mL methanol and shake
- Place centrifuge tube for 10 min in a ultra-sonic bath
- Shake the tube and repeat 5 times
- Centrifuge the mixture at 4500 rpm, for 10 min at 25 °C
- Take 2.5 mL of the centrifugate and dilute it with 2.5 mL of water
- Use the mixture for solid phase extraction

Solid phase extraction

Column:	CHROMABOND [®] PFAS, 6 mL, 300 mg, (REF 730283)
Conditioning:	10 mL 0.1 % NH₃ in methanol, 10 mL methanol, 10 mL water
Sample application:	5 mL of mixture with a flow rate of 2–3 mL/min
Washing:	5 mL of 25 mM ammonium acetate buffer (pH 4.0) with a flow rate of 3 mL/min
Drying:	1 min with vacuum
Elution:	7.5 mL 0.1 % NH_3 in methanol
Eluent exchange:	Evaporate eluate to dryness at 40 °C under a stream of nitrogen and dissolve residue in 0.5 mL water/methanol (20:80, v/v)

Subsequent analysis: HPLC-MS/MS

Chromatographic conditions:

0 1							
Column:	EC 50/2 NUCLEOSHELL [®] RP 18plus, 2.7 µm (REF 763232.20)						
Eluent A:	5 mM ammonium acetate in water						
Eluent B:	5 mM ammonium acetate in methanol						
Gradient:	hold 40 % B for 0.5 min, in 4 min from 40 % B to 95 % B, hold 95 % B for 1.5 min, in 0.05 min to 40 % B, hold 5 % B for 1.45 min						
Flow rate:	0.3 mL/min						
Temperature:	40 °C						

MACHEREY-NAGEL GmbH & Co. KG · Neumann-Neander Str. 6–8 · 52355 Düren · Germany www.mn-net.com · info@mn-net.com

Injection volume:	2 μL	Collision gas:	medium
MS conditions: AB So	ciex QTRAP 5500	lonspray voltage:	-4500 V
Acquisition mode:	SRM	Temperature:	400 °C
Interface:	ESI	lon source gas 1:	50
Polarity:	negative	lon source gas 2:	60
Curtain gas:	30	Detection window:	60 sec

MRM transitions

Abbreviation	Compound	Q ₁	Q ₃ (quan.)	Q₃ (qual.)	RT (min)
3,6-OPFHpA	Perfluoro-3,6-dioxaheptanoic acid	200.9	85.0	134.9	2.18
PFBA	Perfluoro-n-butanoic acid	212.9	168.8	88.9	0.78
PF4OPeA	Perfluoro-4-oxapentanoic acid	228.9	84.9	197.0	0.93
PFPeA	Perfluoro-n-pentanoic acid	262.9	219.0	68.7	1.47
PF50HxA	Perfluoro-5-oxahexanoic acid	279.1	85.0	229.0	1.74
FBSA	Perfluoro-1-butanesulfonamide	297.9	77.9	183.9	2.48
L-PFBS	Perfluoro-1-butanesulfonate	298.9	79.9	98.9	1.64
PFHxA	Perfluoro-n-hexanoic acid	312.9	268.8	119.0	2.29
PFEESA	Perfluoro(2-ethoxyethane)sulfonate	315.1	135.1	69.1	2.00
4:2FTS	1H,1H,2H,2H-Perfluoro-1-hexanesulfonate	326.9	306.9	81.0	2.22
HFPO-DA	2,3,3,3-Tetrafluoro-2-(1,1,2,2,3,3,3-heptafluoropropoxy)propanoic acid	328.9	284.8	169.0	2.49
L-PFPeS	Perfluoro-1-pentanesulfonate	348.9	79.9	98.9	2.36
PFHpA	Perfluoro-n-heptanoic acid	362.9	318.7	169.0	2.86
NaDONA	Sodium dodecafluoro-3H-4,8-dioxanonanoate	376.9	250.7	85.0	2.94
FHEA	2-Perfluorohexyl ethanoic acid	377.0	292.8	95.0	2.96
FHxSA	Perfluoro-1-hexanesulfonamide	398.0	78.0	96.9	3.55
PFHxSK	Perfluoro-1-hexanesulfonate	398.9	79.8	98.9	2.89
PFOA	Perfluoro-n-octanoic acid	412.9	369.0	169.0	3.27
6:2FTS	1H,1H,2H,2H-Perfluoro-1-octanesulfonate	426.9	406.9	79.9	3.25
L-PFHpS	Perfluoro-1-heptanesulfonate	448.9	79.8	98.9	3.28
PFNA	Perfluoro-n-nonanoic acid	462.9	418.9	169.0	3.61
FOEA	2-Perfluorooctyl ethanoic acid	476.9	392.8	412.9	3.72
FOSA	Perfluoro-1-octanesulfonamide	497.9	77.8	63.9	4.17
PFOSK	Perfluorooctanesulfonate	498.8	79.9	99.0	3.60
PFDA	Perfluoro-n-decanoic acid	512.8	468.9	219.1	3.89
8:2FTS	1H,1H,2H,2H-Perfluoro-1-decanesulfonate	526.8	506.8	81.0	3.87
9CI-PF3ONS	9-Chlorohexadecafluoro-3-oxanonane-1-sulfonate	530.8	350.7	82.8	3.78
L-PFNS	Perfluoro-1-nonanesulfonate	548.8	79.9	98.8	3.87
PFUdA	Perfluoro-n-undecanoic acid	562.8	518.9	169.1	4.12
N-MeFOSAA	N-Methylperfluoro-1-octanesulfonamidoacetic acid	569.8	418.9	168.9	4.01
FDEA	2-Perfluorodecyl-ethanoic acid	576.9	493.0	512.8	4.22
N-EtFOSAA	N-Ethylperfluoro-1-octanesulfonamidoacetic acid	583.8	418.8	168.9	4.13
L-PFDS	Perfluoro-1-decanesulfonate	598.8	79.9	98.9	4.09
PFDoA	Perfluoro-n-dodecanoic acid	612.9	568.9	169.0	4.32
11CI-PF3OUdS	11-Chloroeicosafluoro-3-oxaundecane-1-sulfonate	630.7	451.0	82.8	4.25
PFTrDA	Perfluoro-n-tridecanoic acid	662.8	618.9	169.0	4.50
PFTeDA	Perfluoro-n-tridecanoic acid	712.7	668.8	168.9	4.65
6-2diPAP	(1H,1H,2H,2H-Perfluorooctyl)phosphate	788.9	97.0	78.9	4.63
6-2-8-2diPAP	(1H,1H,2H,2H-Perfluorooctyl-1H,1H,2H,2H-perfluorodecyl)phosphate	888.9	442.8	97.1	4.84
8-2diPAP	(1H,1H,2H,2H-Perfluorodecyl)phosphate	988.9	96.9	78.8	5.01
MPFBA	Perfluoro-n-(13C4)butanoic acid	216.9	171.9	-	0.75

Compound	Q ₁	Q ₃ (quan.)	Q_3 (qual.)	RT (min)
Perfluoro-n-(¹³ C ₅)pentanoic acid	267.9	222.9	-	1.47
Perfluoro-1-(2,3,4-13C3)butanesulfonate	301.9	98.9	-	1.64
Perfluoro-n-(1,2-13C2)hexanoic acid	314.9	269.8	-	2.29
1H,1H,2H,2H-Perfluoro-1-(1,2-13C2)hexanesulfonate	328.9	81.0	-	2.21
Perfluoro-n-(1,2,3,4-13C4)heptanoic acid	366.9	321.8	-	2.85
2-Perfluorohexyl-(1,2- ¹³ C ₂)ethanoic acid	379.1	293.8	-	2.96
Perfluoro-1-(1,2,3-13C3)hexanesulfonate	401.9	79.9	-	2.88
Perfluoro-n-(¹³ C ₈)octanoic acid	420.9	376.0	-	3.27
1H,1H,2H,2H-Perfluoro-1-(1,2-13C2)octanesulfonate	428.9	81.0	-	3.25
Perfluoro-n-(¹³ C ₉)nonanoic acid	471.9	427.0	-	3.61
2-Perfluorooctyl-(1,2-13C2)ethanoic acid	478.9	394.0	-	3.73
Perfluoro-1-(13C8) octanesulfonamide	505.9	77.9	-	4.17
Perfluoro-1-(¹³ C ₈)octanesulfonate	506.9	98.9	-	3.60
Perfluoro- <i>n</i> -(1,2,3,4,5,6- ¹³ C ₆)decanoic acid	518.9	474.0	-	3.88
1H,1H,2H,2H-Perfluoro- 1-(1,2-13C2)decanesulfonate	528.9	80.9	-	3.88
Perfluoro- <i>n</i> -(1,2,3,4,5,6- ¹³ C7)undecanoic acid	569.9	525.0	-	4.12
N-Methyl-d3-perfluoro-1-octanesulfonamidoacetic acid	572.9	419.0	-	4.01
2-Perfluorodecyl-(1,2- ¹³ C ₂)ethanoic acid	578.9	493.8	-	4.23
N-Ethyl-d5-perfluoro-1-octanesulfonamidoacetic acid	588.8	418.8	-	4.14
Perfluoro-n-(1,2- ¹³ C ₂)dodecanoic acid	614.9	569.8	-	4.33
Perfluoro-n-(1,2-13C2)tridecanoic acid	714.9	670.0	-	4.66
	CompoundPerfluoro- n -(${}^{13}C_{5}$)pentanoic acidPerfluoro- n -($1,2,3,4-{}^{13}C_{3}$)butanesulfonatePerfluoro- n -($1,2,{}^{13}C_{2}$)hexanoic acid1H,1H,2H,2H-Perfluoro- $1-(1,2,{}^{13}C_{2})$ hexanesulfonatePerfluoro- n -($1,2,3,4-{}^{13}C_{4}$)heptanoic acid2-Perfluorohexyl-($1,2-{}^{13}C_{2}$)ethanoic acidPerfluoro- n -($(1,2,3,{}^{-13}C_{4})$ heptanoic acidPerfluoro- n -($(1,2,{}^{-13}C_{2})$ ethanoic acidPerfluoro- n -($({}^{13}C_{8})$ octanoic acidPerfluoro- n -(${}^{13}C_{9}$)octanoic acid2-Perfluorooctyl-($(1,2-{}^{13}C_{2})$ ethanoic acidPerfluoro- n -(${}^{13}C_{8}$)octanoic acid2-Perfluorooctyl-($(1,2-{}^{13}C_{2})$ ethanoic acidPerfluoro- n -(${}^{13}C_{8}$)octanesulfonatiePerfluoro- n -(${}^{13}C_{8}$)octanesulfonatiePerfluoro- n -($1,2,3,4,5,6-{}^{13}C_{6}$)decanoic acid1H,1H,2H,2H-Perfluoro- $1-(1,2-{}^{13}C_{2})$ decanesulfonatePerfluoro- n -($1,2,3,4,5,6-{}^{13}C_{7}$)undecanoic acid N -Methyl- d_{3} -perfluoro- $1-$ octanesulfonamidoacetic acid2-Perfluorodecyl-($1,2-{}^{13}C_{2}$)ethanoic acid N -Ethyl- d_{5} -perfluoro- $1-$ octanesulfonamidoacetic acid2-Perfluoro- $n-(1,2,{}^{13}C_{2})$ ethanoic acid N -Ethyl- d_{5} -perfluoro- $1-$ octanesulfonamidoacetic acidPerfluoro- $n-(1,2-{}^{13}C_{2})$ dodecanoic acidPerfluoro- $n-(1,2-{}^{13}C_{2})$ tridecanoic acidPerfluoro- $n-(1,2-{}^{13}C_{2})$ tridecanoic acid	CompoundQ1Perfluoro- n -(${}^{13}C_{8}$)pentanoic acid267.9Perfluoro-1-(2,3,4- ${}^{13}C_{3}$)butanesulfonate301.9Perfluoro- n -(1,2- ${}^{13}C_{2}$)hexanoic acid314.91H,1H,2H,2H-Perfluoro-1-(1,2- ${}^{13}C_{2}$)hexanesulfonate328.9Perfluoro- n -(1,2,3,4- ${}^{13}C_{4}$)heptanoic acid366.92-Perfluoro-hexyl-(1,2- ${}^{13}C_{2}$)ethanoic acid379.1Perfluoro-1-(1,2,3- ${}^{13}C_{3}$)hexanesulfonate401.9Perfluoro- n -(${}^{13}C_{9}$)octanoic acid420.91H,1H,2H,2H-Perfluoro-1-(1,2- ${}^{13}C_{2}$)octanesulfonate428.9Perfluoro- n -(${}^{13}C_{9}$)nonanoic acid471.92-Perfluoro-1-(${}^{13}C_{9}$)octanesulfonatic505.9Perfluoro-1-(${}^{13}C_{9}$)octanesulfonatic506.9Perfluoro-1-(${}^{13}C_{9}$)octanesulfonate506.9Perfluoro- n -(1,2,3,4,5,6- ${}^{13}C_{9}$)decanoic acid518.91H,1H,2H,2H-Perfluoro-1-(1,2- ${}^{13}C_{9}$)decanoic acid572.92-Perfluoro- n -(1,2,3,4,5,6- ${}^{13}C_{7}$)undecanoic acid572.92-Perfluoro- n -(1,2,3,4,5,6- ${}^{13}C_{7}$)undecanoic acid578.9N-Methyl-d_{3}-perfluoro-1-octanesulfonamidoacetic acid578.9N-Ethyl-d_{5}-perfluoro-1-octanesulfonamidoacetic acid588.8Perfluoro- n -(1,2,- ${}^{13}C_{2}$)docdecanoic acid614.9Perfluoro- n -(1,2,- ${}^{13}C_{2}$)docdecanoic acid714.9	CompoundQ, <th< td=""><td>Compound Q1 Q3 (qual.) Q3 (qual.) Perfluoro-n-(¹³C₉)pentanoic acid 267.9 222.9 - Perfluoro-n-(¹³C₉)pentanoic acid 301.9 98.9 - Perfluoro-n-(1,2¹³C₉)pexanoic acid 314.9 269.8 - 1H, 1H,2H.2H-Perfluoro-1-(1,2¹³C₉)pexanesulfonate 328.9 81.0 - Perfluoro-n-(1/2,3,4¹³C₉)pexanesulfonate 366.9 321.8 - 2-Perfluorohexyl-(1,2¹³C₉)pexanesulfonate 401.9 79.9 - Perfluoro-n-(¹³C₉)pexanesulfonate 401.9 79.9 - Perfluoro-n-(¹³C₉)pexanesulfonate 420.9 376.0 - 1H,1H,2H,2H-Perfluoro-1-(1,2¹³C₉)pexanesulfonate 428.9 81.0 - Perfluoro-n-(¹³C₉)pexanesulfonate 428.9 81.0 - Perfluoro-n-(¹³C₉)pexanoic acid 471.9 427.0 - 2-Perfluoro-n-(¹⁴C₉)catanesulfonamide 505.9 77.9 - Perfluoro-n-(¹⁴C₉)catanesulfonate 506.9 98.9 - Perfluoro-n-(1,²,3,4,5,6-¹³C₉)decanoic a</td></th<>	Compound Q1 Q3 (qual.) Q3 (qual.) Perfluoro-n-(¹³ C ₉)pentanoic acid 267.9 222.9 - Perfluoro-n-(¹³ C ₉)pentanoic acid 301.9 98.9 - Perfluoro-n-(1,2 ¹³ C ₉)pexanoic acid 314.9 269.8 - 1H, 1H,2H.2H-Perfluoro-1-(1,2 ¹³ C ₉)pexanesulfonate 328.9 81.0 - Perfluoro-n-(1/2,3,4 ¹³ C ₉)pexanesulfonate 366.9 321.8 - 2-Perfluorohexyl-(1,2 ¹³ C ₉)pexanesulfonate 401.9 79.9 - Perfluoro-n-(¹³ C ₉)pexanesulfonate 401.9 79.9 - Perfluoro-n-(¹³ C ₉)pexanesulfonate 420.9 376.0 - 1H,1H,2H,2H-Perfluoro-1-(1,2 ¹³ C ₉)pexanesulfonate 428.9 81.0 - Perfluoro-n-(¹³ C ₉)pexanesulfonate 428.9 81.0 - Perfluoro-n-(¹³ C ₉)pexanoic acid 471.9 427.0 - 2-Perfluoro-n-(¹⁴ C ₉)catanesulfonamide 505.9 77.9 - Perfluoro-n-(¹⁴ C ₉)catanesulfonate 506.9 98.9 - Perfluoro-n-(1, ² ,3,4,5,6- ¹³ C ₉)decanoic a

Table 1: MRM transitions and retention times of PFAS.



Figure 2: Chromatogram of PFAS standard solution after eluent exchange procedure on NUCLEOSHELL® RP 18plus column (β = 2.5 ng/mL for each compound).



Figure 3: Chromatogram of PFAS extract from sample matrix sand on NUCLEOSHELL[®] RP 18plus column (β = 2.5 ng/mL for each compound).



Figure 4: Chromatogram of PFAS from sample matrix soil on NUCLEOSHELL® RP 18plus column (β = 2.5 ng/mL for each compound).

Calibration and Recovery rates

Abbreviation	% Reco- very (sand)	% RSD	% Reco- very (so <u>il)</u>	% RSD	R ₂	Abbreviation	% Reco- very (sand)	% RSD	% Reco- very (so <u>i</u> l)	% RSD	R ₂
11CI-PF3OUdS	84.6	3.5	64.2	2.4	0.995	PFHxA	84.6	3.7	75.2	4.2	0.9
3.6-OPFHpA	92.9	3.2	81.0	6.8	0.993	PFHxSK	90.5	2.3	74.4	3.6	0.9
4:2FTS	91.3	4.7	71.2	5.1	0.995	PFNA	85.1	4.3	58.4	8.8	0.9
6:2FTS	133.9	33.7	58.5	9.3	0.996	PFOA	73.5	3.4	48.2	1.9	0.9
6-2-8-2diPAP	74.1	8.2	29.5	7.8	0.994	PFOSK	86.3	3.0	74.4	7.9	0.9
6-2diPAP	84.4	8.8	39.5	8.1	0.996	PFPeA	90.7	3.2	80.6	6.3	0.9
8:2FTS	90.4	3.9	58.1	12.5	0.997	PFTeDA	89.8	3.4	45.0	6.2	0.9
8-2diPAP	64.2	11.1	29.3	10.1	0.993	PFTrDA	84.8	3.0	49.0	6.1	0.9
9CI-PF3ONS	86.2	5.1	66.6	4.1	0.997	PFUdA	84.0	3.2	55.4	5.1	0.9
FBSA	38.9	15.4	29.7	52.7	0.992	d3- <i>N</i> -MeFO- SAA	66.7	7.3	29.0	6.5	0.9
FDEA	46.2	6.5	34.3	21.3	0.995	d5-N-EtFOSAA	70.3	15.0	35.1	9.9	0.9
FHEA	29.9	12.5	42.8	7.3	0.991	M2-4FTS	89.0	5.3	72.5	5.2	0.9
FHxSA	20.7	11.8	28.9	38.3	0.993	M2-6:2FTS	99.6	7.6	60.8	9.3	0.9
FOEA	27.4	12.4	30.6	6.6	0.994	M2-8:2FTS	97.3	5.8	53.2	8.8	0.9
FOSA	46.2	10.8	31.9	27.5	0.996	M2PFTeDA	89.2	4.4	47.0	4.8	0.9
HFPO-DA	85.0	4.5	82.0	7.4	0.994	M3PFBS	92.3	3.9	84.5	2.9	0.9
L-PFBS	93.9	3.4	85.7	5.3	0.993	M3PFHxS	90.3	3.8	73.6	5.8	0.9
L-PFNS	87.5	1.5	62.3	3.7	0.997	M4PFBA	89.5	2.8	85.3	3.0	0.9
L-PFPeS	90.0	3.7	80.6	2.2	0.994	M4PFHpA	88.7	3.7	70.4	1.9	0.9
NaDONA	85.3	4.0	70.1	3.5	0.995	M5PFHxA	90.6	2.9	77.4	3.9	0.9
N-EtFOSAA	73.0	6.9	33.1	7.7	0.996	M5PFPeA	92.8	3.7	81.5	2.0	0.9
N-MeFOSAA	69.3	4.7	30.1	10.9	0.994	M6PFDA	90.3	4.2	56.1	6.0	0.9
PF4OPeA	89.2	4.3	79.8	2.3	0.993	M7PFUdA	88.3	2.6	55.1	5.4	0.9
PF50HxA	87.1	7.5	79.9	3.9	0.993	M8FOSA	42.5	11.1	30.1	26.6	0.9
PFBA	88.5	3.8	82.2	2.6	0.994	M8PFOA	86.9	3.1	66.4	3.0	0.9
PFDA	87.1	5.2	55.7	4.5	0.996	M8PFOS	82.9	5.0	61.4	11.4	0.9
PFDoA	89.8	2.6	52.8	4.5	0.995	M9PFNA	88.4	6.1	60.5	2.6	0.9
L-PFDS	81.2	6.1	59.6	5.5	0.995	MFDEA	44.9	13.4	31.4	9.7	0.9
PFEESA	93.2	2.8	84.6	1.9	0.995	MFHEA	31.5	8.8	39.6	17.3	0.9
PFHpA	90.2	2.2	71.1	3.7	0.993	MFOEA	30.1	12.3	34.3	8.3	0.9
L-PFHpS	85.5	3.9	72.5	7.0	0.995	MPFDoA	89.3	4.1	53.1	4.3	0.9

Table 2: Recovery rates for the presented SPE method for contaminated soils. Correlation coefficient is given for all compound from calibration curves, determined of 7 levels between 0.1 ng/mL and 10.0 ng/mL for each compound.





Figure 5: Recovery rates of PFAS from contaminated soils (spiked with 5 ng/g sample for each compound) according to DIN 3814-14.



Figure 6: Recovery rates of 40 PFAS from contaminated soils (spiked with 5 ng/g sample for each compound).

Conclusion

This application note shows the reliable and successful determination of per- and polyfluoroalkyl substances (PFAS) from sediments like sand and soil with an optimized SPE method. The successful determination of 40 PFAS was successfully carried out with the special CHROMABOND[®] PFAS column. In this way, the combination of different SPE sorbents in a multi-layer column allows to use various of interaction types like ionic, hydrophobic, hydrogen bonds, dipole-dipole and π - π interactions for the enrichment of a broad spectrum of PFAS in a single enrichment procedure.

With the presented methodology good recovery rates for PFAS from matrices like sand and sufficient recovery rates from soil with good reproducibility could be achieved. Figure 5 shows that high recovery rates for PFAS analysis according to DIN could be obtained. More results for the additional 30 PFAS are presented in table 2 and figure 6.

The investigation of the extraction procedure described in DIN shows that the efficiency of the extraction of PFAS depends on different parameters like polarity and adsorption effects. These effects lead to losses of PFAS before the SPE enrichment procedure. Probably a second extraction step with a mixture of methanol and water could lead to better extraction efficiency of more PFAS.

The chromatographic separation of PFAS was performed by using core-shell particles that are well known for fast and high-efficient separations combined with a reasonably low backpressure. In this work, subsequent analysis was developed on a NUCLEOSHELL[®] RP 18plus column as shown in figure 2. The chromatographic results provide a good correlation for all PFAS compounds as presented in table 2.

In summary, the presented application describes a quick and convenient method for the determination of various PFAS from sediment samples with a SPE procedure using the extraction mechanisms of DIN 30407-42.

References

- Solid phase extraction of per- and polyfluoroalkyl substances (PFAS) from drinking water, MN application note 05/2020
- [2] German standard methods for the examination of water, waste water and sludge - Jointly determinable substances (group F) – Part 42: Determination of selected polyfluorinated compounds (PFC) in water – Method using high performance liquid chromatography and mass spectrometric detection (HPLC/MS-MS) after solid-liquid extraction (F 42), 2011–03.

Product information

The following MACHEREY-NAGEL products have been used in this application note:

REF 763232.20	EC	50/2	HPLC	column	(analytical),
	NUCL	EOSHEL	L [®] RP 18p	olus, 2.7 µm	I
REF 730283	CHRO	DMABON	D [®] PFAS,	6 mL, 300 i	mg
REF 702402	Screv Polyin	v closure, nide oranç	N 9, PP, ge, 1.0 mn	blue, c. h n, flourine-f	ole, Sili.w./ ree
REF 702009	Screv inner	v neck via cone, PP	al, N 9, 1 tr.	1.6x32.0 m	nm, 0.3 mL,

MACHEREY-NAGEL shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance, or use of this material. Information, descriptions, and specifications in this publication are subject to change without notice.