Size Exclusion Columns for Antibodies




YMC-SEC MAB

Features

- Especially designed for antibody analysis
- Best choice for Bio QC
- Simultaneous analysis of monomers, fragments and aggregates
- Excellent resolution and peak shapes suitable


Excellent Lot-to-Lot Reproducibility



10 μL (5 mg/mL) humanised monoclonal antibody Injection:


YMC-SEC MAB provides excellent reproducibility of the separation of monomer and aggregates as well as for monomer and fragments. Therefore, it is very effective for quality control of antibody drugs.

Ideal Choice for Monoclonal Antibodies




Simultaneous analysis of both high and low molecular weight regions is possible with just one column: YMC-SEC MAB.

Suitable for Antibody-Drug-Conjugates (ADCs)



YMC-SEC MAB is also suitable for the analysis of Antibody-Drug Conjugates (ADC). The addition of an organic solvent to the mobile phase is effective in improving the results obtained for ADC analysis.

Ideal for Analysis of Digested Antibodies



Column: YMC-SEC MAB (3 μ m, 25 nm) 300 x 4.6 mm ID Part No.: DLM25S03-3046WT


Eluent: 0.1 M phosphate buffer (pH = 7) cont. 0.2 M NaCl Flow rate: 0.165 mL/min

Flow rate: 0.165 mL/min
Temperature: 25 °C
Detection: UV at 280 nm
Injection: 2 µL (3 mg/mL)

 $\begin{array}{ll} \mbox{Injection:} & 2 \ \mu\mbox{L (3 mg/mL)} \\ \mbox{Sample:} & \mbox{Humanised monoclonal IgG1 + Papain} \end{array}$

Ideal for the analysis of fragments/degradation products of antibodies:
Digestion of a monoclonal antibody with papain was monitored for 7 hours. The peak of the monomer decreased as digestion proceeded, while peaks for degradation products increased.

High Column Stability



Excellent stability is provided for monoclonal antibody analysis without any changes in theoretical plate number or elution time even after more than 500 injections.

Specifications			
MW range	10 – 700 kDa		
Particle sizes	3 μm (<2 μm*)		
Pore size	25 nm / 250 Å		
Max. Temperature	40 °C		
Max. Pressure	14 MPa (2,030 psi)		
Base particle	silica / dihydroxypropyl		

^{*}scheduled for a later release

Ordering Information

Particle size (µm)	Pore size (nm)	Column ID (mm)	Column Length (mm)	Part number
3	25	4.6	150	DLM25S03-1546WT
		4.6	300	DLM25S03-3046WT
		8.0	300	DLM25S03-3008WT