



**NUCLEOSIL® SA**  
**NUCLEOSIL® SB**

**Note:** All HPLC columns from MACHEREY-NAGEL are supplied with a certificate, which contains specifications and test results of the column. The columns NUCLEOSIL® SA and NUCLEOSIL® SB are quality products based on the robust silica NUCLEOSIL®. They are specifically developed for HPLC analysis. If carefully and properly used excellent chromatographic results and long column lifetime can be achieved. HPLC columns are designed for qualitative and quantitative analysis of mixtures of substances and single components. They must exclusively be used in accordance with universally accepted laboratory regulations and HPLC working methods. Before running the column the entire analytical system (column and equipment) has to be carefully checked by the operator. Chromatographic conditions (mobile phase, flow, temperature etc.) must be adapted to the analytical task. MACHEREY-NAGEL does not give any warranty and is not liable for the success of a separation or application. If you have any questions after reading this leaflet, please call our service / technical support.

**Table of contents**

- Safety indication
- Description of the column
- Installation
- Guard columns
- Sample
- Eluent
- Flow rate and pressure
- Temperature
- Detection
- Equilibration
- Column storage
- Troubleshooting
- Column regeneration
- Abstract

**Safety indication**

Follow the general safety instructions for handling of HPLC solvents used as mobile phases (e.g., additives of methanol, acetonitrile or acids) and take precautions against any kind of injuries or damage to health (e.g., skin and eye protection in case of broken capillaries). Disposal of used HPLC columns must follow international, national and local environmental protection regulations. The use of HPLC columns is only permitted to staff members, who are qualified in their field. Keep HPLC columns away from children. MACHEREY-NAGEL disclaims and excludes all warranties of any kind or nature whatsoever and MN shall not be liable for any damages (whether direct, indirect, foreseeable, incidental, compensatory, consequential or special), whether based upon warranty, contract, tort or strict liability, if damages and/or losses occur caused by improper use, maintenance, neglect or improper treatment (especially opening of the column and exposure of the column bed).

**Description of the column**

As stationary phases NUCLEOSIL® SA columns contain a strongly acidic cation exchanger (SCX), made of spherical silica modified with benzenesulfonic acid. NUCLEOSIL® SB columns are packed with a strongly basic anion exchanger (SAX), a quaternary ammonium modification. Both phases have an ion exchange capacity of ~ 1 meq/g.

Typical applications of the cation exchanger NUCLEOSIL® SA are separations of basic, water soluble substances. The retention of basic analytes is governed by pH value and ion strength of the mobile phase (see eluent). The anion exchanger NUCLEOSIL® SB is used for the separation of substances, which form anions in aqueous solutions. Such substances are mainly organic acids (aromatic and aliphatic carboxylic and sulfonic acids) as well as pesticides and pharmaceuticals.

**Installation**

The columns should be installed in the flow direction indicated on the column label. They are connected with 1/16" capillaries and fittings, typical for HPLC instruments.

**Guard columns**

For protection and an extension of column lifetime columns should always be used with a guard column. The filter elements and the adsorbent in the guard column retain contaminants from the sample or the eluent. Connection of the guard column with the separation column is made by a suitable guard column holder (see [www.mn-net.com](http://www.mn-net.com) or the MN chromatography catalog). Cartridge replacement is required when increased column pressure and/or loss of performance is observed.

**Sample**

Sample solutions should be passed through a syringe filter (e.g., CHROMAFIL® Xtra PET, 0.45 µm, 25 mm, REF 729220) before entering the column. If injected sample solutions are still turbid even after filtration, the lifetime of the column may be significantly reduced. The sample volume should be as small as possible to achieve an optimal resolution.

**Eluent**

The ion exchanger columns NUCLEOSIL® SA and SB are supplied with a solution of 0.15 mol/L diammonium hydrogen phosphate, pH 5. As mobile phase aqueous buffers (e.g., phosphate, acetate, citrate) are used. Eluents should be filtered through a 0.2–0.45 µm membrane filter and degassed. The pH value can be adjusted in relation to the pK<sub>s</sub> value of the analyte with the correspondent acid (e.g., phosphoric, acetic, citric acid). The pH should be in the range of 2–8. Strongly acidic or basic conditions can result in dissolution of the column bed or the organic modification. The amount of buffer salts should be kept as low as possible (max. buffer concentration: 0.15 mol/L). Organic additives (e.g., methanol, acetonitrile) can be used to improve sample solubility or the separation. Take care of a possible precipitation of buffer salts and plugging of the column when adding an organic additive.

**Flow rate and pressure**

Flow rate (recommended for analytical columns with 2–4.6 mm ID: 0.2–2.0 mL/min) influences the time required, the resolution and the column lifetime. It is limited by the back pressure, which should not exceed the maximum of 400 bar. We recommend controlling back pressure regularly. If a high pressure results from the use of the column at nominal flow rates, this usually indicates that some contaminants have become deposited on the packing material, which must be removed (see troubleshooting).

**Temperature**

Column temperatures up to 60 °C are possible. For a long lifetime, especially under basic conditions, 30–40 °C is recommended. However, they should be at least 30 °C below the boiling temperature of the eluent, in order to ensure proper detection. Variation of the temperature influences retention times and especially the peak shape. Optimum temperatures for successful separations should be determined empirically.

**Detection**

Spectrophotometers, refractometers and electrochemical detectors can be used with the columns. If electrochemical detectors are used, please note that high temperatures may be incompatible with some working electrodes. If a higher sensitivity is required, post-column derivatizations with an appropriate detector for the reaction product can be used.

**Equilibration**

Prior to measurement of samples the column must be rinsed with the eluent at the same flow rate and temperature as the method to be applied. Column equilibration is finished, when the baseline of the detector no longer shows a drift (generally after 10 column volumes).

**Column storage**

The original eluent 0.15 mol/L diammonium hydrogen phosphate, pH 5 is recommended for storage. To avoid bacterial growth, store the column in a refrigerator. For column storage be sure the end fittings are tightly sealed using column end plugs, because storage without these seals can result in drying of the packing material. Under these circumstances rinse the column with approx. 10 column volumes of the eluent of storage at a flow rate of max. 0.2 mL/min.

**Troubleshooting**

The following outline describes the symptoms of performance loss and their cause. All columns are subject to the strict regulation and control of our quality assurance system. Columns based on silica are robust and hold their separation efficiency for long periods by correct maintenance and treatment. According to experience, column failures are mostly a result of injection of contaminants to the sorbent bed. The usage of a guard column, as well as an appropriate sample pretreatment will help to minimize these risks.

Use the outline below to help determine the cause of a possible performance loss:

Symptom / Error / Cause	Prevention / Remedy
<b>Baseline drift</b> • insufficient period for equilibration with the eluent • contaminated eluent • temperature	longer or better equilibration use freshly prepared solvents and reagents column temperature control
<b>Broad peaks</b> • mixing and/or diffusion before/behind the column • too large sample volume	keep length and ID of capillaries at a minimum smaller injection volume
<b>Peak interference; too fast elution</b> too fast elution and/or insufficient separation by: • improper column temperature or flow rate • elution power of eluent is too high	optimize concerned parameter optimize eluent system
<b>Increasing back pressure; degradation of the separation performance</b> contamination of sorbent by: • particulate accumulation on frit or sorbent bed from sample, eluent or system • precipitation of buffer salts	prepare fresh eluent; prefilter samples and eluent, use in-line filter / rinse LC system, clean the sorbent check solubility of buffer salts before / remove them by rinsing (see column regeneration)
<b>Insufficient separation; degradation of the separation with regular column pressure</b> contamination with: • fats, oils, lipids from sample (coating of sorbent surface) and other organic substances from improperly prepared eluent or matrices	remove organic substances by sample preparation / clean the sorbent (see column regeneration)
<b>Double peaks (dead volume)</b> • faulty fittings (capillaries, ferrules, nuts) • dissolution of silica by too high pH value of eluent	use "PEEK Fingertight Fittings", REF 718770 / replace fittings consider pH range of column / replace column

**Column regeneration**

In some cases the performance of the column can be restored by removing contaminants from the sorbent bed or by regeneration of the phase. It is important, however, to locate the source of contamination before using the column for the analysis of samples again.

- Prepare fresh eluent:** In some cases the performance loss is traced to eluent contamination. Therefore, prepare fresh eluent and flush all liquid lines before using the column again. The eluent should be filtered through a 0.2–0.45 µm membrane and degassed prior to use.
- Cleaning of sorbent:** To remove contamination rinse the column with a minimum of 10 column volumes (see table below) at the original flow rate and temperature as follows:
  - 100% water for removal of the buffer and of ionic impurities (Cationic impurities on the SA phase can be also removed with 1 mol/L sodium perchlorate, pH <4)
  - 100% methanol to remove nonionic organic compounds (remove buffer before)
  - if necessary, 100% methanol with inverse flow direction at 1/5 of original flow rate
  - rinse the column with 10 column volumes 100% water at original flow direction
  - convert column to storage condition with 0.15 mol/L diammonium hydrogen phosphate, pH 5 at original flow rate

An adequate indicator for a clean column is a constant baseline. At constant temperature you should observe less than 2–3 mAU drift during a running time of 5 minutes with an isocratic run.

- Column replacement:** The above procedures will restore performance only in certain cases. Some organic contaminants are particularly refractory and may not respond to treatment. Also dead volume, due to column compression can generally not be repaired. Under these circumstances, column replacement is necessary. It is highly advisable to locate the cause of the problem before installing a new column.



Length [mm]	Inner diameter [mm]:	Column volume [mL]			
		2	3	4	4.6
100		0.30	0.70	1.25	1.65
150		0.45	1.05	1.90	2.50
250		0.80	1.75	3.15	4.15

**Abstract**

To extend column lifetime, please keep in mind the following:

- As eluents aqueous buffers (e.g., phosphate, acetate, citrate) are recommended. Avoid a precipitation of buffer salts when using of organic additives (e.g., methanol). Eluents should be filtered through a 0.2–0.45 µm membrane and degassed.
- Filter samples through a 0.2–0.45 µm CHROMAFIL® Xtra PET syringe filter before injection.
- Use a guard column for contaminated samples.
- The recommended flow rate for analytical columns (ID 2–4.6 mm) is 0.2–2.0 mL/min.
- Adjust flow rate to keep column pressure below 400 bar.
- Store the column in 0.15 mol/L diammonium hydrogen phosphate, pH 5 in a refrigerator.
- Use analytical grade reagents and HPLC grade solvents for all work. Discard any solutions that show evidence of bacterial growth.

**Please check the full range of MACHEREY-NAGEL chromatography products!**



... for applicative support please ask for our HPLC Application Guide or visit our website with more than 3000 chromatography applications: [www.mn-net.com/apps](http://www.mn-net.com/apps)