

6.4 Volt Electronic Crimper & Decapper



Features:

- Powerful motor with a fast cycle
- Long lasting battery with Lithium Ion technology
- Advanced charging system including warning light when recharging is necessary
- Can be operated while charging
- Convenient controls on top of unit to adjust crimp
- Optional base and mounting kit (5ABAS0)



CRS electronic crimpers and decappers fit standard laboratory vials and seals with aluminum or two-part aluminum/ steel caps. The electronically adjustable cycle is accurate, reproducible, and can be adjusted with the control pad on top of the unit.

The tools are ergonomically designed to reduce strain and arm injury associated with many manual crimping tools. Available in 8 (crimper only), 11, 13, and 20 mm cap sizes.

Rechargeable

The tools can be recharged in one to two hours with the universal power supply included. The power supply includes a plug set for operation in most countries.



Description	Item No.
8 mm Electronic Crimper	5A08C0
11 mm Electronic Crimper	5A11C0
13 mm Electronic Crimper	5A13C0
20 mm Electronic Crimper	5A20C0
13 mm Electronic Crimper, for Flip Off Caps	5A13F0
20 mm Electronic Crimper, for Flip Off Caps	5A20F0

Description	Item No.
11 mm Electronic Decapper	5A11D0
13 mm Electronic Decapper	5A13D0
20 mm Electronic Decapper	5A20D0
Base with Mounting Kit	5ABAS0

