

# Manual Hand Crimper & Decapper




## Innovative Design Body

The body is constructed of tough, light-weight glass-filled plastic for strength. CRS's ergonomically curved handle design conforms to your natural hand grip, which greatly improves comfort over blocky metal grips.

The bottom-pull handle design is steadier to hold than top-push handle designs. The head does not wobble. It is easy to hold the tapered jaws of the crimper in a stable position, as the bottom handle is pulled up.

## Adjustment Knob

- Easily-viewed
- Easy to see when crimp setting has been reached



## Tough Glass-filled Plastic Body

- Strong light-weight material

## Bottom-Pull Handle Design


- Stable position design compared to top-push handle designs

## Curved High-Strength Plastic Handles

- Conforms to your natural grip
- Improving comfort over blocky metal grips

## Adjusting the Crimp Setting

The adjustment knob is conveniently positioned at the head of the crimper, clearly displays + and - symbols along with directional arrows. This simplifies adjusting the intensity of crimp desired. No more "extra squeeze" necessary! Allows to clearly see the crimp setting has been reached when the knob touches the crimper body.



To ensure a smooth and consistent operation, the hand crimper uses the same precision tolerances as the CRS Electronic Crimper and Decapper.

Description	Item No.
11 mm Easy Grip Manual Crimper	311990
20 mm Easy Grip Manual Crimper	320990
13 mm Easy Grip Manual Crimper	313990
8 mm Easy Grip Manual Crimper	308990

Description	Item No.
13 mm Easy Grip Manual Crimper, for Flip Off Caps	313992
20 mm Easy Grip Manual Crimper, for Flip Off Caps	320992
11 mm Easy Grip Manual Decapper	311991
20 mm Easy Grip Manual Decapper	320991
13 mm Easy Grip Manual Decapper	313991

