



La-Pha-Pack Stainless Steel Cleanroom Crimping Tools

Ensure sample integrity in cleanroom environments

La-Pha-Pack cleanroom crimping tools are made entirely of stainless steel for maximum durability during dry heat sterilization or autoclaving.

- The crimping/decapping mechanism is corrosion and heat resistant while the stainless steel construction removes the need for any protective coating on the handle or crimp head.
- Can withstand repeated sterilization for cleanroom use without the risk of damaging the tool.
- Available in 11mm, 13mm and 20mm sizes.
- Crimping tools are adjustable in crimping pressure and height to offer optimal crimping results on varying vial styles.

Ordering information:

La-Pha-Pack Cleanroom Crimping Tools

Description	Size	Part No
Stainless steel crimper	11mm	11 06 2225
Stainless steel crimper	13mm	13 06 2227
Stainless steel crimper	20mm	20 06 2229
Stainless steel decapper	11mm	11 07 2226
Stainless steel decapper	13mm	13 07 2228
Stainless steel decapper	20mm	20 07 2230



Other sizes are available upon request.

For more information, please contact customer service:

Tel.: +49 2423 94 31 – 0

Fax: +49 2423 94 31 – 34

e-mail: service@la-pha-pack.com

www.la-pha-pack.com

©2012 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

Part of Thermo Fisher Scientific