



La-Pha-Pack

OEM Solutions for Chromatography



Your reliable partner
for Chromatography Consumables

Total Phthalat Free (TPF) Seals with assembled aluminum liner

- Septa material free of any elastomers and halogens; therefore no contamination of the sample by plasticizers, Silicone or Butyl rubbers and PTFE / FEP / TEF materials
- Recommended for GC/MS and LC/MS applications
- The additional ring which is assembled above the aluminum liner achieves an excellent tightness and inertness of the complete combination seal
- Suitable for storing of standards or reactive substances which would normally corrode the septum material

Application areas:

- Analytic of elastomers und plastomers
- Phthalate analysis
- Analytic of fluorinated, halogenated organic compounds
- VOC (volatile organic compounds) analysis
- Analytic of polymerisation catalysts
- Analytic of silicones und silicates

Ordering information:

La-Pha-Pack TPF Seals

Description	Article No.
11mm Combination Seal: Aluminum Cap, clear lacquered, centre hole, aluminum liner, 0.06mm (Sealing by an additional installed transparent silicone-flat seal)	11 03 2519
9mm Combination Seal: PP Short Thread Cap, blue, centre hole; aluminium liner, 0.06mm (Sealing by an additional installed transparent silicone-flat seal)	09 15 2487



La-Pha-Pack GmbH
Am Parir 20 – D-52379 Langerwehe
Phone: +49 (0) 24 23 / 9431-0 – Fax: +49 (0) 24 23 / 9431-33
Email: service@la-pha-pack.com – Internet: www.la-pha-pack.com

ThermoFisher
SCIENTIFIC